MS 10310030114

Desproteinizador

Eletrólitos

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ENGLISH

IMPORTANCE OF USE

The use of the deproteinizing solution maintains the performance of the K, Ca and Cl selective ion electrodes of the ELETROLYTE - WE-300 instrument from WAMA Diagnóstica.

Its use is necessary for a high performance of the K, Ca and Cl selective ion electrodes whenever the instrument alerts for the electrode maintenance.

The **Electrolyte Deproteinizer** from **WAMA** uses the pepsin enzyme and a diluter containing potassium chloride, sodium chloride and hydrochloric acid, thereby activating the response of the K, Ca and Cl electrodes of the ELETROLYTE - WE-300 when this becomes slow or unstable.

PRINCIPLE OF THE METHOD

The ELETROLYTE - WE-300 uses the selective ion electrode principle to determine the ion concentrations of Na^+ , K^+ , Cl^- , Ca^{2+} , in addition to pH, in serum, plasma, and urine. The proteins and fat present in the serum form a barrier on the surface of the K^+ , Ca^{2+} and Cl^- electrodes, diminishing its performance after a long period of use. The deproteinizing solution removes this material and reactivates the performance of the electrodes.

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REF 836005-B – 5 ampoules

- 1- Diluter: 5 x 0,8ml
- 2- Enzyme: 5 x 50mg
- 3- Instructions for use

STORAGE AND REAGENT STABILITY

1 - Electrolyte Deproteinizer: is supplied in 2 vials, one diluter (0.8 mL) and one enzyme (50 mg), and must be kept in cool, dry and away from intense light, and the temperature must be between 2-8 °C. Do not freeze. The expiration date of the reagents is printed on the label of the vials. Once open, it must be used immediately.

PROCEDURE

1. Add the contents of the ampoule of the diluent (0.8 mL) into the

ampoule of the enzyme. Shake the bottle until the solution is transparent.

2. The solution is stable for 20 to 30 minutes, so it must be used before this time.

3. Follow the instrument's User's Manual ELETROLYTE WE-300 for aspiration or transferring the solution and continuity of the procedure.

PRECAUTIONS AND WARNINGS

1. For in vitro Diagnostic use only.
2. Read carefully the instructions for use prior to performing the test.
3. The expiry date refers to the last day of the month indicated on the labels of the vial and of the kit box.
4. Avoid exposing the reagents to high temperatures and direct sunlight.
5. Do not freeze the reagent, as this may cause irreversible damage.
6. Avoid contact with skin and eyes.
7. Do not use kit components after the expiration date.
8. Discard the material following local regulations.
9. Use the Good Laboratory Practices (GLPs) in storage, handling and disposal of materials.

WARRANTY

WAMA Diagnóstica guarantees the exchange of the diagnostic kit, provided it is within the expiration date and is proven by its technical support that there were no technical mistakes in the execution, handling and storage of this product. WAMA and its distributors are not liable for failures on the performance of the kits under these conditions.

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BIBLIOGRAFIA / BIBLIOGRAPHY / BIBLIOGRAFIA

1. Heins, M., Heil, W. and Withold, W.: Storage of serum or whole blood samples. Effects of time and temperature on 22 serum analytes, Eur. J. Clin. Chem. Clin. Biochem., (33): 221 – 228, 1995.
2. Burnett, R.W. et al.: Approved IFCC recommendations on whole blood sampling, transport and storage for simultaneous determination of pH, blood gases and electrolytes., Eur. J. Clin. Chem. Clin. Biochem., (33): 247 – 253, 1995.
3. Young, D.S.: Effects of drugs on clinical laboratory tests. 5th Ed. Washington DC: AACC Press; 2000.
4. Young, D.S.: Effects of preanalytical variables on clinical laboratory tests. 3rd Ed. Washington DC: AACC Press; 2007.

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